

Review article

Current and Future Percutaneous Strategies for the Treatment of Acute and Chronic Heart Failure

Ignacio J. Amat-Santos,^{a,*} Luis H. Varela-Falcón,^a and William T. Abraham^b^a Instituto de Ciencias del Corazón (ICICOR), CIBERCV, Hospital Clínico Universitario de Valladolid, Valladolid, Spain^b Cardiology Department, Ohio State University, Columbus, Ohio, United States

Article history:

Available online 30 January 2017

Keywords:

Heart failure
Cardiogenic shock
Percutaneous devices

Palabras clave:

Insuficiencia cardiaca
Shock cardiogénico
Dispositivos percutáneos

ABSTRACT

The prevalence of heart failure (HF) has risen in parallel with improved survival in patients after a myocardial infarction and an aging population worldwide. In recent years, new percutaneous therapies have been developed to complement current established treatments for acute/decompensated and chronic HF and minimize risks. In acute presentations, the failure of medical treatment is no longer the end of the road in refractory circulatory shock; the use of mechanical circulatory support devices may be the next milestone in well-resourced health settings. Although evidence in this area is difficult to generate, research networks can facilitate the volume and quality of data needed to further augment the clinician's knowledge. Pulsatile (intra-aortic balloon pump), axial continuous (Impella), or centrifugal continuous pumps (TandemHeart; HeartMate PHP) together with percutaneously implanted extracorporeal membrane oxygenation are radically changing the prognosis of acute HF. Newer percutaneous therapies for chronic HF are based on attractive hypotheses, including left atrial decompression with shunting devices, left ventricle restoration through partitioning devices, or pressure-guided implantable therapies that may help to promptly treat decompensations. To date, only the last has been proved effective in a randomized study. Therefore, thorough research is still needed in this dynamic and promising field.

© 2016 Sociedad Española de Cardiología. Published by Elsevier España, S.L.U. All rights reserved.

Terapias percutáneas en el tratamiento de la insuficiencia cardiaca aguda y crónica: presente y futuro

RESUMEN

La prevalencia de la insuficiencia cardiaca (IC) aumenta paralelamente al envejecimiento poblacional y la mayor supervivencia de los pacientes con infarto de miocardio. Recientemente se han desarrollado nuevas terapias percutáneas con el objetivo de complementar los tratamientos actuales de la IC aguda/descompensada y crónica minimizando los riesgos. En la IC aguda, el fracaso del tratamiento médico ya no puede ser el final de la estrategia terapéutica para el *shock* cardiogénico, dado el éxito de los dispositivos de soporte circulatorio mecánico. A pesar de que la evidencia en esta área es difícil de generar, el pronóstico de la IC aguda está cambiando radicalmente gracias a la investigación en red sobre dispositivos de flujo pulsátil (balón de contrapulsación intraaórtico), continuo axial (Impella) o continuo centrífugo (TandemHeart; HeartMate PHP) junto con el implante percutáneo de oxigenador extracorpóreo de membrana. Las nuevas terapias percutáneas en la IC crónica se basan en atractivas hipótesis como la descompresión de la aurícula izquierda (*shunts*), la restauración ventricular (compartimentación) o la monitorización intratorácica de presiones mediante dispositivos implantables que permitan tratar precozmente las descompensaciones. Actualmente solo la última se ha probado efectiva en un estudio aleatorizado. Por lo tanto, es necesario estudiar a fondo este dinámico y prometedor campo.

© 2016 Sociedad Española de Cardiología. Publicado por Elsevier España, S.L.U. Todos los derechos reservados.

* Corresponding author: Instituto de Ciencias del Corazón (ICICOR), Hospital Clínico Universitario de Valladolid, Ramón y Cajal 3, 47005 Valladolid, Spain.
E-mail address: ijamat@gmail.com (I.J. Amat-Santos).

Abbreviations

AMI: acute myocardial infarction
 HF: heart failure
 LA: left atrium
 LV: left ventricular
 LVAD: left ventricular assist device
 MCS: mechanical circulatory support
 NYHA: New York Heart Association

INTRODUCTION

The prevalence of heart failure (HF) has risen to 2% in adults in parallel with its progressive increase during the life course.¹ World population aging has led to a growing global incidence of HF conditions, impairing patient wellbeing and increasing health care costs.^{1,2} Several improvements in the treatment of patients with HF with reduced ejection fraction have reduced the mortality rate, which continues to be higher than that of acute myocardial infarction (AMI), and the condition is associated with a rehospitalization rate of 40% per year.^{1,2} Hence, the mid-term prognosis of HF is still dire and in many scenarios is similar to that of cancer.² In addition, HF represents between 1.8% and 3.1% of the public health budget in our environment, mostly due to in-hospital costs (73% in Spain).²

Therefore, efforts should focus on reducing the incidence of HF and decompensations, to avoid hospitalizations, and on improving prognosis. Several pharmacological and nonpharmacological strategies, including defibrillators, cardiac resynchronization therapy, left ventricular (LV) assist devices (LVAD), and heart transplant, represent important milestones in the treatment of advanced HF.^{1,3,4} In recent years, new percutaneous therapies have been developed to simplify and reduce the aggressiveness of interventions to complement current established treatments for acute/decompensated and chronic HF. In this review, we aim to summarize the main current and developing percutaneous alternatives for the treatment of both acute and chronic HF.

PERCUTANEOUS THERAPIES FOR ACUTE HEART FAILURE

Cardiogenic shock is defined as a systemic tissular hypoperfusion secondary to low cardiac output despite adequate circulatory volume and LV filling pressure. As a consequence, mean blood pressure persists below 90 mmHg or decreases by 30 mmHg or more, with a cardiac index below 1.8 L/min/m² (without hemodynamic support) or below 2.2 L/min/m² if the patient is under support, and a wedge pulmonary pressure \geq 15 mmHg.^{5–8}

The mortality rate in this serious scenario has improved in the last few decades but continues to be unacceptably high (\approx 50%), as does morbidity from the resultant multiple organ failure. The failure of inotropic and vasoactive drugs to overcome this critical status is no longer the end of the road in refractory circulatory shock; the use of percutaneous mechanical circulatory support (MCS) devices may be the next milestone in well-resourced health settings. The aim of MCS is two-fold: to stabilize the patient and to evaluate the subsequent options that can range from bridge to decision, to recovery, to long-term support devices such as percutaneous ventricular assist devices or to total artificial hearts, and/or heart transplant.⁵ The choice of the initial rescue MCS strategy has a significant bearing not only on limiting further

iatrogenic harm in the acute setting, but also on planning long-term strategies in the absence of myocardial recovery. The development of predetermined institutional pathways is critical to the success of such MCS programs. Indeed, the ideal timing for percutaneous MCS device implantation remains controversial and the only general consensus is that they should be implanted sooner rather than later. In advanced HF, the INTERMACS score, which includes 7 categories, helps to better identify those stages (1 and 2) that may benefit from MCS due to failed inotropics therapy.^{9,10} Expert consensus is that percutaneous ventricular assist devices are indicated in patients undergoing high-risk coronary angioplasty (class IIB) and in those with cardiogenic shock due to AMI (class IB),^{5,11} although their potential use in cardiogenic shock due to other etiologies is also recognized (Table 1 of the supplementary material).

Although evidence in this area is difficult to generate, research networks can facilitate the volume and quality of data needed to further augment the clinician's knowledge of when, where, and which type of these technologies could and should be used. Mechanical circulatory support devices are classified according to the type of blood pump: Pulsatile (intra-aortic balloon pump [IABP]), axial continuous (Impella), or centrifugal continuous (TandemHeart; CentriMag; Rotaflow; HeartMate PHP).¹² The IABP remains the most widely used method for mechanical assistance⁵ in patients experiencing LV failure, because of its moderate hemodynamic capabilities, prompt time to therapy, and relatively low complication rates. Percutaneous LVADs (pLVADs) represent an emerging option for partial or total circulatory support and several studies have compared the safety and efficacy of these devices with those of IABP. A description of each technique is provided below and summarized in Table 1.






Intra-aortic Balloon Pump (Arrow IABP, Winston-Salem, North Carolina, United States)

Intra-aortic balloon pump is still the most widely used percutaneous ventricular assist device worldwide even though there is no proven evidence of its clinical benefit.¹³ Nevertheless, American and European guidelines recommend its use in patients with AMI and shock (class IIA),^{23,24} and in 2014 half of the patient receiving a heart transplant in Spain were under hemodynamic support with IABP.^{25,26} The main components include a 7.5 to 8.0 Fr dual-lumen catheter and a polystyrene balloon inflated with helium that allows rapid inflation and deflation due to its low viscosity. This creates a pulsatile pumping action with the deflation occurring during systole, which helps to reduce the afterload, increasing coronary perfusion and decreasing oxygen demand by the myocardium, but with a modest increase in cardiac output.²⁷ Therefore, below a certain threshold of cardiac output, it may not be effective. In addition, its use is precluded by electrical instability, vascular disease, or the presence of aortic regurgitation.

Impella (Abiomed, Danvers, Massachusetts, United States)

A few randomized controlled trials have demonstrated better hemodynamic profiles for pLVADs compared with IABP.^{16,14,28} However, this has not translated into improved 30-day survival. Moreover, patients treated with pLVADs tended to have a higher incidence of leg ischemia and device-related bleeding.²⁹ The use of positive inotropic drugs or vasopressors was expected to be lower in patients with Impella but no differences have been detected in the overall use of these agents. However, the use of Impella devices may increase in patients not responding to inotropes and IABP. Therefore, the decision-making process on treatment requires an integrated stepwise approach. In 2012 O'Neill et al.¹⁶ reported the

Table 1
Technical Characteristics and Main Studies of Current Percutaneous Mechanical Circulatory Support Devices

Type of device	IABP (pulsatile flow)	Impella (axial continuous flow)	HeartMate PHP (centrifugal continuous flow)	TandemHeart (centrifugal continuous flow)	VA-ECMO (centrifugal continuous flow)
Figures					
Size, Fr	7-8	2.5 4 (CP) 5	13	V- 21 A- 15-19	V- 18-31 A- 15-23
Output, LPM	0.5	2.5 4 5	4	3.5-4.5	4.5-8
Complexity of implantation	+	++	+++ (surgical)	++	+++
Potential complications	Limb ischemia Bleeding	Limb ischemia Bleeding Hemolysis	Limb ischemia Bleeding Hemolysis	Limb ischemia Bleeding Hemolysis	Limb ischemia Bleeding Hemolysis
Trials	1. SHOCK-II ¹³ (vs medical treatment)	1. ISAR shock ¹⁴ (Impella 2.5 vs IABP) 2. RECOVER right trial ¹⁵ (Impella CP) 3. PROTECT-II trial ¹⁶ (Impella 2.5 vs IABP) 4. IMPRESS trial ¹⁷ (Impella CP vs IABP)	1. SHIELD-I ¹⁸ 2. SHIELD-II ¹⁹	1. TH study ²⁰ (vs IABP) 2. TH registry ²¹	1. Registry of ECMO in cardiac arrest ²² 2. Lack of randomized trials
Clinical settings	1. Cardiogenic shock	1. Cardiogenic shock 2. RV failure after 48 h with LVAD 3. High-risk PCI	1. High-risk PCI 2. High-risk PCI	1. Cardiogenic shock 2. Right heart failure	1. Out-of-hospital cardiac arrest
Results	1. No differences in 30-d mortality (48% IABP; 49% medical)	1. No difference in 30-d mortality (46% both groups) 2. 30-d mortality: 27% 3. No differences in 30-d mortality (40.1% IABP; 35.1% Impella)	1. Safety study (not randomized) 2. Ongoing	1. No differences in 30-d mortality (45% IABP; 43% TH) 2. In-hospital mortality: 57%	1. In-hospital mortality: 50%

A, arterial; ECMO, extracorporeal membrane oxygenation; IABP, intra-aortic balloon pump; LPM, liters per minute; LVAD, left ventricular assist device; PCI, percutaneous coronary intervention; RV, right ventricle; TH, TandemHeart; V, venous; 30-d, 30-day.

Images reproduced with permission from the copyright owner.

results of the PROTECT II study, a prospective, randomized trial comparing hemodynamic support with Impella 2.5 vs IABP in individuals undergoing high-risk percutaneous coronary intervention. The trial was discontinued prematurely due to futility. Currently, there is also limited evidence on whether Impella improves the net health outcome in individuals with cardiogenic shock in the postcardiotomy setting or in high-risk percutaneous coronary interventions. Early study findings suggest that, although hemodynamic measures are consistently greater, clinically meaningful changes in outcomes have not been demonstrated.

Two percutaneous Impella devices exist: the Impella 2.5 and the Impella CP, both with a pigtail shape distally, which is introduced in the left ventricle, where a continuous flow is created toward the ascending aorta.³⁰ Due to this mechanism, the flow does not require an activation signal—unlike the IABP—and is not affected by the presence of arrhythmias. Currently, it is approved for use for up to 4 (Food and Drug Administration [FDA]) or 5 days (*Conformité Européenne* [CE] mark).³¹ As limitations, it requires an adequately functioning right ventricle or the use of a right ventricular assist device, and its use is contraindicated by the presence of mechanical aortic prostheses or severe vascular disease. Complications are also similar to those of IABP, including bleeding events related to the need for systemic anticoagulation, and vascular injury. More specifically, inadvertent malpositioning of the device is a particular concern that may lead to rapid deterioration, and hemolysis may occur in up to 10% of patients.^{30–32} Both issues can be solved, in general, by repositioning and pump impeller adjustment.

TandemHeart (CardiacAssist Inc, Pittsburgh, Pennsylvania, United States)

This device has a 21-Fr cannula allocated in the left atrium (LA) through transeptal puncture impelling blood from the LA to the iliofemoral artery through a 15 to 19-Fr cannula.³⁰ Its use has been approved by the FDA for 6 hours whereas a CE mark has been obtained for up to 30 days.

Like the Impella device, right ventricular function should be preserved³⁰ and vascular disease is a contraindication. There is limited information concerning its use in patients with aortic regurgitation or mechanical complications related to AMI.^{33,34} A specific source of complications is the need for transeptal puncture that should be guided by transesophageal echocardiography to minimize the risk.⁵ In addition, potential displacement of the cannula allocated in the LA to the right atrium is an infrequent but challenging event. Several studies comparing TandemHeart to IABP revealed that the cardiac index rose in both groups (and in general hemodynamics improved), but was significantly higher in the TandemHeart group.^{35,36} In contrast, overall mortality at 30 days was similar in both groups, whereas the rate of events such as leg ischemia, severe bleeding, or sepsis was higher among TandemHeart participants. As before, the small number of patients probably did not allow for a meaningful evaluation of potential mortality differences but the results suggest that TandemHeart should only be used if the patient's hemodynamic status is not sufficiently improved by IABP.

Extracorporeal Membrane Oxygenation

Recently, there has been increasing application of extracorporeal membrane oxygenation (ECMO) technology to provide MCS in an incremental fashion either as peripheral or central venoarterial ECMO or venovenous assist devices. The use of ECMO in cardiopulmonary resuscitation is also expanding with experienced centers reporting favorable outcomes.^{37,38} The percutaneous

approach has been proved to be safe, effective and time saving, while minimizing bleeding from the cannula insertion site. Moreover, ECMO can be implanted at the bedside without the need for fluoroscopy guidance.^{5,30} However, percutaneous cannulation remains a challenge and is not free of complications. Indeed, ischemia of the corresponding extremity still requires surgical revision in 17% of patients. Therefore, a careful insertion technique, close surveillance, and monitoring are essential. The size of the cannula is similar to that of TandemHeart and portable devices are relatively simple to implant. The most frequent complications are related to thrombosis, bleeding,³⁹ and infective events. As with TandemHeart, an antegrade perfusion cannula is recommended to prevent limb ischemia. The FDA and CE mark have approved the continuous use of the ECMO for 30 days, but longer duration has been consistently reported by the registry of the Extracorporeal Life Support Organization.²²

There is a lack of randomized trials but the consensus is that ECMO can be useful in cases of cardiogenic shock, postcardiotomy failure, myocarditis, during resuscitation,^{37,38} or if there is primary graft failure in heart transplants. Indeed, one third of patients who underwent emergent heart transplant in Spain in 2014 were under ECMO.¹⁴

One of the main reasons for increased morbidity and mortality in patients under MCS is right ventricular failure. Currently, the use of venoarterial ECMO is probably the best alternative if this complication occurs, although the results of the RECOVER RIGHT study suggest a new potential and interesting use for Impella in this setting.^{15,38}

Newer Percutaneous Therapies for Acute Heart Failure







The HeartMate PHP (Percutaneous Heart Pump) from Thoratec (St. Jude) is a new continuous axial flow catheter whose insertion is similar to that of the Impella device (Table 1) and with CE mark approval since July 2015 for its use in high-risk percutaneous coronary interventions. Its potential use in cardiogenic shock will be evaluated by the SHIELD-II study (coronary interventionS in High-risk patiEnts using a novel percutaneous LV support device). Theoretically, it may diminish the rate of hemolysis through a reduction in shear stress and increase the stability within the left ventricle.

Several promising new technologies in preclinical research stages are also of interest in this context. This is a highly dynamic field and 2 examples that represent only the tip of the iceberg are the Aortix (Procyon, Houston, Texas, United States) and the BoLetz micro-LVAD (OCR, Yale, United Kingdom). The first consists of a self-expanding frame that fixes the pump to the aortic wall within the descending aorta, avoiding the risk of emboli to the supra-aortic vessels; following implantation, the 18-Fr catheter is removed, leaving only a flexible electrical power wire which is tunneled to a transcutaneous energy transfer system. The endovascularly deployable BoLetz micro-LVAD is a new concept in the early stages of development, which may represent the future concept of pLVAD destination therapies.

PERCUTANEOUS THERAPIES FOR CHRONIC HEART FAILURE

One of every 5 people will suffer from HF during their lifetime.^{1,2} One of the main reasons for this is the current better survival after an AMI, which has led to a higher prevalence of chronic LV dysfunction. However, HF still carries a poor prognosis and, in particular, progressive HF represents the main cause of death in patients in New York Heart Association (NYHA) class III or IV. Therefore, limiting this progression is a major target to improve

Table 2
New Percutaneous Therapies for the Treatment of Chronic Heart Failure

Therapeutic strategies	LA decompression			LV restoration	Pressure-monitoring	Injectable therapies	Others
Devices	AFR (Mia Medical, Istanbul, Turkey)	V-Wave (V-Wave Ltd, Or-Akiva, Israel)	IASD (DC Devices Inc, Tewksbury, Massachusetts, United States)	Parachute (CardioKinetix Inc, Menlo Park, California, United States)	CardioMEMS HF (CardioMEMS, Atlanta, Georgia, United States)	Cell/Gene therapy + Helix transcatheter delivery system (BioCardia, San Carlos, California, United States)	New resynchronization devices, Neuromodulation, etc.
Figures							NA
Size, Fr	14	14	16	14-16	12	Variable	-
Complexity of implantation		++		++++	++	+++	+++
Trials/studies	Registry (PHT)	Registry (HFpEF) ⁴¹	Registry (HFpEF) ⁴²	1. Registries ⁴³ 2. PARACHUTE-IV ⁴⁴ (trial ongoing)	1. Registries ⁴⁵ 2. CHAMPION trial ⁴⁶	Registries	Registries
Results	Reduction in rehospitalizations			LV volume reduction and 6MWT distance improvement	Reduction in rehospitalizations. Included in HF guidelines	Under evaluation	Under evaluation
Other devices	"Home-made" devices			VenTouch (Mardil Medical, Minneapolis, Minnesota, United States)	Chronicle (Medtronic, Minneapolis, Minnesota, United States)	Algisyl (LoneStar Heart Inc, Laguna Hills, California, United States)	Multiple

6MWT, 6-minute walk test; HF, heart failure; HFpEF, heart failure with preserved ejection fraction; HFrEF, heart failure with reduced ejection fraction; LA, left atrial; LV, left ventricular; NA, not available; PHT, pulmonary hypertension.

Images reproduced with permission from the copyright owner.

not only quality of life but also mid-term survival.¹ In addition, although there have been major developments in therapies for chronic HF with reduced ejection fraction, very few advances have been made in the heterogeneous scenario of HF with preserved ejection fraction, which lacks therapies with proven effectiveness, despite its comparable rates of morbidity and mortality.¹ Given the large number of patients affected by this problem and the huge costs for the health system, it is not surprising that prolific research has been performed in this field in recent decades.⁴⁰ The main new percutaneous strategies for the treatment of chronic HF are schematically depicted in Table 2. Percutaneous valvular therapies, as well as those aimed to treat mechanical complications, are not the object of this review.

Most HF patients have symptoms related to congestion and to low cardiac output. New therapies focus mainly on reducing left chamber filling pressures to improve congestion and decrease myocardial oxygen consumption.^{45,47} With this aim, different mechanisms have been proposed, including pressure-monitoring implantable devices, left atrial decompression, and LV reshaping systems.

Pressure-monitoring Devices

For decades now, we have been using direct intracardiac pressure measurements through a Swan-Ganz catheter to obtain accurate diagnosis, evaluate prognosis, and treat patients with HF, based on the findings of classic retrospective studies. However, prospective studies including SUPPORT,⁴⁸ ESCAPE,⁴⁹ or PAC-Man,⁵⁰ have failed to demonstrate its benefits but invasive pressure monitoring is still a commonly used tool in the management of HF patients.^{51,52} The missing gap between clinical practice and current evidence may be met by implantable pressure-monitoring devices. This new technology aims to prevent and promptly treat decompensations in chronic HF, through intracardiac wireless devices that can register LV filling pressures.⁵³ The Chronicle device (Medtronic, Minneapolis, Minnesota, United States) makes a continuous recording of this pressure through a wire similar to pacemaker leads, within the right ventricular outflow tract. The COMPASS-HF study⁵⁴ demonstrated that clinical decompensation was preceded by an increase in filling pressures, and helped to achieve a 36% reduction in the relative risk of rehospitalization in patients with chronic HF in NYHA class III. Similar findings were obtained from the REDUCE HR study,⁵⁵ which was, however, discontinued prematurely due to a wire malfunctioning. For this reason, the FDA refused to approve its clinical use in 2007 but it nevertheless represents the proof of concept for the strategy of wireless pressure monitoring.

More recently, the CHAMPION study,^{45,46} which included HF patients in NYHA class III randomly assigned to conventional or pressure-guided treatment, reported that the CardioMEMS HF device (CardioMEMS, Atlanta, Georgia, United States) achieved a 28% reduction in the rehospitalization rate within the first 6 months following implantation. This device has been approved by the regulatory authorities, leading to its inclusion in the latest clinical guidelines for patients with HF.¹ This small monitor (15 × 3 mm) has a silicon encapsulation and nitinol anchoring to avoid its migration. Implantation is performed through a 12-Fr catheter advanced to the pulmonary artery that can be charged through radiofrequency.

Finally, the HeartPOD (St. Jude Medical) monitors LA pressure through a 3 × 7 mm monitor anchored in the interatrial septum. The ongoing LAPTOT-HF and CRT-D LAP studies⁵⁶ aim to evaluate LA pressures and subsequently inform the patient on how to titrate the medications accordingly. The main limitation of this new device may be the need for interatrial septal puncture and its associated risks but the risk/benefit ratio is still under evaluation.

Left Atrial Decompression Shunting Devices

The rationale for shunting blood between atria as a means to reduce elevated pressures is based on medical knowledge from the 1960s (creation of atrial septostomies) to divert blood from one atrium to the other in the treatment of elevated pressures (both right and left off-loading).⁵⁷ Moreover, currently the use of septostomy procedures is recommended by pulmonary hypertension guidelines for patients with advanced pulmonary hypertension with elevated right atrial pressure, and specific devices such as the Heartx AFR have been designed to this end (Table 2). Published data and feasibility studies in HF animal models demonstrate that shunting of 1000 to 1500 mL/min for a 6 mmHg pressure gradient (between left and right arterial pressure) and 1500 to 1900 mL/min for a 10 mmHg pressure gradient is expected to reduce LAP by 20% and 30%, respectively, without significantly increasing right atrial pressure.^{58,59} The population for which these new shunting devices are indicated includes patients with clinical sequels due to high LA pressure, representing approximately 90% of the patients hospitalized for pulmonary congestion.⁴⁷ The reduction of this pressure below 18 mmHg could potentially improve the patient's lung congestive state with a direct and positive effect on ventilation, dyspnea, and hospitalization rates.^{41,42,60,61}

Two small series have reported the use of specific left-to-right shunting devices in humans to date. The first represents initial experience with the V-Wave device (V-Wave Ltd, Or Akiva, Israel) for the treatment of HF with reduced ejection fraction.⁴¹ The V-Wave device is a dedicated device, designed for percutaneous implantation in the fossa ovalis. This device attempts to address some of the limitations characterizing atrial septostomies. It is composed of a self-expandable hour-glass shaped nitinol shunt, with an expanded polytetrafluoroethylene encapsulation and 3 porcine pericardial leaflets sutured together. The V-Wave delivery catheter is a 14-Fr catheter intended to deliver the crimped shunt. Oral anticoagulation is recommended for 3 months afterward. The device was successfully implanted in 10 patients; no device-related or procedural adverse events occurred during follow-up. The NYHA class improved from III to II-I in 89% of the patients, with a significant increase in quality of life according to specific questionnaires and in the 6-minute walk test distance. Pulmonary capillary wedge pressure was reduced from 23 ± 5 mmHg to 17 ± 8 mmHg at 3 months, with no changes in right atrial pressure, pulmonary arterial pressure, or pulmonary resistance. No patient was admitted to hospital for worsening HF.⁴¹

In the second series, 68 patients with a diagnosis of HF with preserved ejection fraction were treated with the IASD (Corvia Medical, Tewksbury, Massachusetts, United States) with successful placement in 64 patients.⁴² The IASD is also a dedicated device composed of nitinol with an outer and inner diameter of 19 mm and 8 mm, respectively. Thus, a permanent 8-mm atrial septal defect is created. The legs of the device are flat on the LA side to minimize the risk of thrombus formation. A 16-Fr sheath is placed in the femoral vein and a proprietary delivery catheter is used to deliver the implant to the desired location at the fossa ovalis. After the procedure, patients were treated with dual antiplatelet therapy but no anticoagulation was recommended as there is no valve sutured inside the nitinol structure. No patient had periprocedural cardiac or cerebrovascular events or required a cardiac surgical intervention for device-related complications during 6 months of follow-up. A total of 52% had a reduction in pulmonary capillary wedge pressure at rest and 58% during exertion, with an increase in the mean exercise duration, which suggests that atrial blood shunting could be a new strategy for the management of HF with preserved ejection fraction.⁴²

To date, the effectiveness of LA decompression shunting devices has not been compared with the standard of care in a randomized

controlled trial. Therefore, this intriguing hypothesis and preliminary data on hemodynamic and functional improvement need to be interpreted with caution.

Left Ventricular Restoration Percutaneous Devices

The current aggressive invasive strategies for AMI have increased survival rates. However, this has created a growing population with large residual myocardial injury. Left ventricular dilation and remodeling postmyocardial infarction is a well-known phenomenon driven by the reduction in contractility that increases wall stress and promotes the molecular activation of the fibrotic process leading to adverse remodeling. The concept of percutaneous ventricular restoration of the left ventricle, is based on the premise that a dedicated partitioning compliant device delivered via a catheter-based approach may achieve LV volume reduction, geometric reconfiguration, and synchronized wall motion to achieve a more effective ejection while minimizing the risk of a surgical approach. Its therapeutic efficacy depends on improvements in LV volumes reducing the wall stress.

The Parachute system (CardioKinetix Inc, Menlo Park, California, United States) includes a delivery system with a balloon that facilitates expansion of the device, and a preshaped delivery catheter. The device, with alternative sizes, comprises a self-expanding nitinol frame, an expanded polytetrafluoroethylene impermeable membrane with an atraumatic radiopaque polymer foot aimed to provide a contact point between the LV apical wall, which allows the device to be positioned with a vector toward the outflow tract.

After initial promising preclinical results,⁶² a prospective nonrandomized study (PARACHUTE III) aimed to assess the long-term safety and efficacy of the Parachute device. One hundred patients with NYHA class II–IV ischemic HF and LV ejection fraction between 15% and 40% with dilated akinetic or dyskinetic anterior-apical wall were enrolled. Of the 100 enrolled patients, device implantation was successful in 97 participants. The primary safety endpoint was procedural-related major adverse cardiac/cerebral events (7%) and the secondary safety endpoint was the composite of mortality and morbidity (32.3%). Secondary efficacy endpoints included several imaging and functional parameters and were successful in terms of LV volume reduction and improvement in 6-minute walk distance.⁴³ This acceptable safety and efficacy result also needs to be validated in randomized controlled trials, such as the ongoing PARACHUTE IV.⁴⁴

In summary, these novel therapies aiming to improve the symptoms and prognosis of patients with chronic HF are promising, but there is still a gap in the evidence that will soon be narrowed by ongoing studies. The first step, demonstrating the importance of pressure-guided treatment of chronic HF through randomized analysis has already been taken and represents the first milestone for this new generation of percutaneous therapies.

FINAL COMMENTS

Heart failure has a growing incidence and variable mechanisms that require specific therapies. In acute or decompensated presentations, the failure of medical treatment no longer means that therapeutic options have been exhausted in refractory circulatory shock; the use of MCS devices represents the next milestone in well-resourced health settings. Although evidence on survival is scarce because this particular scenario hinders randomization, the growing range of alternative devices must be handled with caution but also bravely in order to reduce the unacceptable mortality rate in patients with cardiogenic shock.

Even more complex is the management of patients with chronic HF. Because of the large number of individuals with this problem, the deterioration in their quality of life, and the health care-related costs, new and equally effective and efficient strategies are needed to improve the quality of life of our patients with a parallel reduction in hospitalizations and decompensations. The slow process of evaluation of new percutaneous therapies to treat chronic HF is both necessary and sufficient to integrate those technologies with demonstrated benefits in our therapeutic arsenal.

In conclusion, the armamentarium of new percutaneous devices for the treatment of acute and chronic HF is leading to a paradigm shift from a common therapeutic path to an individualized approach in the management of these patients. This dynamic field still requires thorough investigation and caution, but also open-minded physicians ready to move forward into this thrilling new era.

ACKNOWLEDGEMENTS

To Amada Recio-Platero for her kind help and to the rest of the professionals of the Heart Failure Unit in the ICICOR (*Instituto de Ciencias del Corazón*).

FUNDING

L.H. Varela-Falcón received support from the Fundación Carolina-BBVA. The ICICOR received support to a related project from the *Instituto de Salud Carlos III* (FIS Project PI 15/01695) and by ERDF funds.

CONFLICTS OF INTEREST

I.J. Amat-Santos has received institutional grant from V-Wave Ltd. W.T. Abraham is consultant and shareholder of V-Wave Ltd and has been a consultant to and received grants from Medtronic Inc, CardioMEMS, St. Jude Medical, and CVRx.

SUPPLEMENTARY MATERIAL



Supplementary material associated with this article can be found in the online version available at [doi:10.1016/j.rec.2016.12.008](https://doi.org/10.1016/j.rec.2016.12.008).

REFERENCES

1. Ponikowski P, Voors AA, Anker SD, et al. ESC guidelines for the diagnosis and treatment of acute and chronic heart failure: The Task Force for the diagnosis and treatment of acute and chronic heart failure of the European Society of Cardiology (ESC). *Eur Heart J*. 2016. <http://dx.doi.org/10.1093/eurheartj/ehw128>.
2. Sayago-Silva I, Gacía-López F, Segovia-Cubero J. Epidemiología de la insuficiencia cardiaca en España en los últimos 20 años. *Rev Esp Cardiol*. 2013;66:649–656.
3. Anguita Sánchez M, Lambert Rodríguez JL, Bover Freire R, et al. Classification and Quality Standards of Heart Failure Units: Scientific Consensus of the Spanish Society of Cardiology. *Rev Esp Cardiol*. 2016;69:940–950.
4. Comin-Colet J, Enjuanes C, Lupón J, Cainzos-Achirica M, et al. Transitions of Care Between Acute and Chronic Heart Failure: Critical Steps in the Design of a Multidisciplinary Care Model for the Prevention of Rehospitalization. *Rev Esp Cardiol*. 2016;69:951–961.
5. Rihal CS, Naidu SS, Givertz MM, et al. 2015 SCAI/ACC/HFSA/STS Clinical Expert Consensus Statement on the Use of Percutaneous Mechanical Circulatory Support Devices in Cardiovascular Care (Endorsed by the American Heart Association, the Cardiological Society of India, and Sociedad Latino Americana de Cardiología Intervención; Affirmation of Value by the Canadian Association of Interventional Cardiology–Association Canadienne de Cardiologie d'intervention). *Catheter Cardiovasc Interv*. 2015;85:175–196.

6. Antonelli M, Levy M, Andrews PJ, et al. Hemodynamic monitoring in shock and implications for management. International Consensus Conference, Paris, France, 27–28 April 2006. *Intensive Care Med.* 2007;33:575–590.
7. Reynolds HR, Hochman JS. Cardiogenic shock: Current concepts and improving outcomes. *Circulation.* 2008;117:686–697.
8. Hasdai D, Topol EJ, Califf RM, et al. Cardiogenic shock complicating acute coronary syndromes. *Lancet.* 2000;356:749–756.
9. Kirklin JK, Naftel DC, Pagani FD, et al. Seventh INTERMACS annual report: 15,000 patients and counting. *J Heart Lung Transplant.* 2015;34:1495–1504.
10. Drakos SG, Kfoury AG, Selzman CH, et al. Left ventricular assist device unloading effects on myocardial structure and function: Current status of the field and call for action. *Curr Opin Cardiol.* 2011;26:245–255.
11. Levine GN, Bates ER, Blankenship JC, et al. 2011 ACCF/AHA/SCAI Guideline for percutaneous Coronary Intervention. A report of the American College of Cardiology Foundation/American Heart Association Task Force on Practice Guidelines and the Society for Cardiovascular Angiography and Interventions. *J Am Coll Cardiol.* 2011;58:44–122.
12. Moazami N, Fukamachi K, Kobayashi M, et al. Axial and centrifugal continuous-flow rotary pumps: A translation from pump mechanics to clinical practice. *J Heart Lung Transplant.* 2013;32:1–11.
13. Thiele H, Zeymer U, Neumann FJ, et al. Intraaortic balloon support for myocardial infarction with cardiogenic shock. *N Engl J Med.* 2012;367:1287–1296.
14. Seyfarth M, Sibbing D, Bauer I, et al. A randomized clinical trial to evaluate the safety and efficacy of a percutaneous left ventricular assist device versus intra-aortic balloon pumping for treatment of cardiogenic shock caused by myocardial infarction. *J Am Coll Cardiol.* 2008;52:1584–1588.
15. Anderson MB, Goldstein J, Milano C, et al. Benefits of a novel percutaneous ventricular assist device for right heart failure: The prospective RECOVER RIGHT study of the Impella RP device. *J Heart Lung Transplant.* 2015;34:1549–1560.
16. O'Neill WW, Kleiman NS, Moses J, et al. A prospective, randomized clinical trial of hemodynamic support with Impella 2.5 versus intra-aortic balloon pump in patients undergoing high-risk percutaneous coronary intervention: The PROTECT II study. *Circulation.* 2012;126:1717–1727.
17. Ouweeneel DM, Eriksen E, Sjaauw KD, et al. Impella CP Versus Intra-Aortic Balloon Pump in Acute Myocardial Infarction Complicated by Cardiogenic Shock: The IMPRESS trial. *J Am Coll Cardiol.* 2016. <http://dx.doi.org/10.1016/j.jacc.2016.10.022>.
18. Weigelt JC. SHIELD I Trial Results Confirm St. Jude Medical HeartMate PHP Cardiac Assist Device Provides Consistent, Reliable Cardiac Support During High-Risk PCI Procedures. 2015 Oct [accessed 2016 Nov 22]. <http://media.sjm.com/newsroom/news-releases/news-releases-details/2015/SHIELD-I-Trial-Results-Confirm-St-Jude-Medical-HeartMate-PHP-Cardiac-Assist-Device-Provides-Consistent-Reliable-Cardiac-Support-During-High-Risk-PCI-Procedures-0/default.aspx>.
19. Supporting Patients Undergoing High-Risk PCI Using a High-Flow Percutaneous Left Ventricular Support Device (SHIELD II) [accessed 2016 Nov 22]. <http://clinicaltrials.gov/show/NCT02468778>.
20. Thiele H, Sick P, Boudriot E, et al. Randomized comparison of intra-aortic balloon support with a percutaneous left ventricular assist device in patients with revascularized acute myocardial infarction complicated by cardiogenic shock. *Eur Heart J.* 2005;26:1276–1283.
21. TandemHeart Experiences and Methods (THEME Registry) [accessed 22 Nov 2016]. <https://clinicaltrials.gov/ct2/show/NCT02326402>.
22. Paden ML, Conrad SA, Rycus PT, et al. ELSO Registry. Extracorporeal Life Support Organization Registry Report 2012. *ASAIO J.* 2013;59:202–210.
23. Roffi M, Patrono C, Collet JP, et al. 2015 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation Task Force for the Management of Acute Coronary Syndromes in Patients Presenting without Persistent ST-Segment Elevation of the European Society of Cardiology (ESC). *Eur Heart J.* 2016;37:267–315.
24. O'Gara PT, Kushner FG, Ascheim DD, et al. 2013 ACCF/AHA guideline for the management of ST-elevation myocardial infarction: A report of the American College of Cardiology Foundation/American Heart Association Task Force on Practice Guidelines. *J Am Coll Cardiol.* 2013;61:78–140.
25. H'Doubler Jr PB, H'Doubler WZ, Bien RC, et al. A novel technique for intraaortic balloon pump placement via the left axillary artery in patients awaiting cardiac transplantation. *Cardiovasc Surg.* 2000;8:463–465.
26. González-Vilchez F, Segovia Cubero J, Almenar L, et al. Registro Español de Trasplante Cardíaco. XXVI Informe Oficial de la Sección de Insuficiencia Cardíaca y Trasplante Cardíaco de la Sociedad Española de Cardiología (1984–2014). *Rev Esp Cardiol.* 2015;68:1008–1021.
27. Papaioannou TG, Stefanadis C. Basic principles of the intraaortic balloon pump and mechanisms affecting its performance. *ASAIO J.* 2005;51:296–300.
28. Ostadal P, Mlecek M, Holy F, et al. Direct comparison of percutaneous circulatory support systems in specific hemodynamic conditions in a porcine model. *Circ Arrhythm Electrophysiol.* 2012;5:1202–1206.
29. Tayal R, Barvalia M, Rana Z, et al. Totally percutaneous insertion and removal of Impella device using axillary artery in the setting of advanced peripheral artery disease. *J Invasive Cardiol.* 2016;28:374–380.
30. Basra SS, Loyalka P, Kar B. Current status of percutaneous ventricular assist devices for cardiogenic shock. *Curr Opin Cardiol.* 2011;26:548–554.
31. Burzotta F, Carlo Trani C, Doshi SN, et al. Impella ventricular support in clinical practice: Collaborative viewpoint from a European expert user group. *Int J Cardiol.* 2015;201:684–691.
32. Lauten A, Engstrom AE, Jung C, et al. Percutaneous left-ventricular support with the Impella-2.5-assist device in acute cardiogenic shock: Results of the Impella-EUROSHOCK-registry. *Circ Heart Fail.* 2013;6:23–30.
33. Pham DT, Al-Quthami A, Kapur NK. Percutaneous left ventricular support in cardiogenic shock and severe aortic regurgitation. *Catheter Cardiovasc Interv.* 2013;81:399–401.
34. Gregoric ID, Bieniarz MC, Arora H, et al. Percutaneous ventricular assist device support in a patient with a postinfarction ventricular septal defect. *Tex Heart Inst J.* 2008;35:46–49.
35. Kar B, Gregoric ID, Basra SS, Idelchik GM, Loyalka P. The percutaneous ventricular assist device in severe refractory cardiogenic shock. *J Am Coll Cardiol.* 2011;57:688–696.
36. Burkhoff D, Cohen H, Brunckhorst C, et al. A randomized multicenter clinical study to evaluate the safety and efficacy of the TandemHeart percutaneous ventricular assist device versus conventional therapy with intraaortic balloon pumping for treatment of cardiogenic shock. *Am Heart J.* 2006;152:469.e1–e8.
37. Romeo F, Acconcia MC, Sergi D, et al. Percutaneous assist devices in acute myocardial infarction with cardiogenic shock: Review, meta-analysis. *World J Cardiol.* 2016;8:98–111.
38. Ortega-Deballon I, Hornby L, Shemie SD, et al. Extracorporeal resuscitation for refractory out-of-hospital cardiac arrest in adults: A systematic review of international practices and outcomes. *Resuscitation.* 2016;101:12–20.
39. Pieri M, Agracheva N, Bonaveglio E, et al. Bivalirudin versus heparin as an anticoagulant during extracorporeal membrane oxygenation: A case-control study. *J Cardiothorac Vasc Anesth.* 2013;27:30–34.
40. Chen J, Normand SL, Wang Y, et al. National and regional trends in heart failure hospitalization and mortality rates for Medicare beneficiaries, 1998–2008. *JAMA.* 2011;306:1669–1678.
41. Del Trigo M, Bergeron S, Bernier M, et al. Unidirectional left-to-right interatrial shunting for treatment of patients with heart failure with reduced ejection fraction: a safety and proof-of-principle cohort study. *Lancet.* 2016;387:1290–1297.
42. Hasenfuß G, Hayward C, Burkhoff D, et al. A transcatheter intracardiac shunt device for heart failure with preserved ejection fraction (REDUCE LAP-HF): a multicentre, open-label, single-arm, phase 1 trial. *Lancet.* 2016;387:1298–1304.
43. Thomas M, Nienaber CA, Ince H, et al. Percutaneous ventricular restoration therapy using the Parachute device in 100 subjects with ischaemic dilated heart failure: one-year primary endpoint results of PARACHUTE III, a European trial. *EuroIntervention.* 2015;11:710–717.
44. Costa MA, Pencina M, Nikolic S, et al. The PARACHUTE IV trial design and rationale: percutaneous ventricular restoration using the parachute device in patients with ischemic heart failure and dilated left ventricles. *Am Heart J.* 2013;165:531–536.
45. Abraham WT, Adamson PB, Bourge RC, et al. Wireless pulmonary artery haemodynamic monitoring in chronic heart failure: a randomised controlled trial. *Lancet.* 2011;377:658–666.
46. Adamson PB, Abraham WT, Aaron M, et al. CHAMPION trial rationale and design: the long-term safety and clinical efficacy of a wireless pulmonary artery pressure monitoring system. *J Card Fail.* 2011;17:3–10.
47. Ritzema J, Troughton R, Melton I, et al. Physician-directed patient self-management of left atrial pressure in advanced chronic heart failure. *Circulation.* 2010;121:1086–1095.
48. Connors Jr AF, Speroff T, Dawson NV, et al. The effectiveness of right heart catheterization in the initial care of critically ill patients. SUPPORT Investigators. *J Am Med Assoc.* 1996;276:889–897.
49. Binanay C, Califf RM, Hasselblad V, et al. Evaluation study of congestive heart failure and pulmonary artery catheterization effectiveness: the ESCAPE trial. *JAMA.* 2005;294:1625–1633.
50. Harvey S, Harrison DA, Singer M, et al. Assessment of the clinical effectiveness of pulmonary artery catheters in management of patients in intensive care (PAC-Man): a randomised controlled trial. *Lancet.* 2005;366:472–477.
51. Stevenson LW, Tillisch JH. Maintenance of cardiac output with normal filling pressures in patients with dilated heart failure. *Circulation.* 1986;74:1303–1308.
52. Steimle AE, Stevenson LW, Chelmsky-Fallick C, et al. Sustained hemodynamic efficacy of therapy tailored to reduce filling pressures in survivors with advanced heart failure. *Circulation.* 1997;96:1165–1172.
53. Castro PF, Concepcion R, Bourge RC, et al. A wireless pressure sensor for monitoring pulmonary artery pressure in advanced heart failure: initial experience. *J Heart Lung Transplant.* 2007;26:85–88.
54. Bourge RC, Abraham WT, Adamson PB, et al. Randomized controlled trial of an implantable continuous hemodynamic monitor in patients with advanced heart failure: the COMPASS-HF study. *J Am Coll Cardiol.* 2008;51:1073–1079.
55. Adamson PB, Gold MR, Bennett T, et al. Continuous hemodynamic monitoring in patients with mild to moderate heart failure: results of The Reducing Decompensation Events Utilizing Intracardiac Pressures in Patients With Chronic Heart Failure (REDUCEhf) trial. *Congest Heart Fail.* 2011;17:248–254.
56. Maurer MS, Adamson PB, Costanzo MR, et al. Rationale and design of the Left Atrial Pressure Monitoring to Optimize Heart Failure Therapy study (LAPTOP-HF). *J Card Fail.* 2015;21:479–488.
57. Amat-Santos IJ, Rodés Cabau J, López J. Left-to-right interatrial shunt percutaneously implanted devices: a new avenue in the treatment of heart failure. *Interv Cardiol.* 2015;7:443–450.

58. Kaye D, Shah SJ, Borlaug BA, et al. Effects of an interatrial shunt on rest and exercise hemodynamics: results of a computer simulation in heart failure. *J Cardiac Fail.* 2014;20:212–221.
59. Keren G, Nitzan Y, Abraham W, et al. An interatrial shunt device reduces left atrial pressure, improves haemodynamics and prevents adverse remodeling in an ovine heart failure model. *EuroIntervention.* 2013;164:A4.
60. Amat-Santos IJ, Bergeron S, Bernier S, et al. Left atrial decompression through left-to-right interatrial shunt for the treatment of left cardiac heart failure: first-in-man experience with the V-Wave device. *EuroIntervention.* 2015;10:1127–1131.
61. SØndergaard L, Reddy V, Kaye D, et al. Transcatheter treatment of heart failure with preserved or mildly reduced ejection fraction using a novel interatrial implant to lower left atrial pressure. *Eur J Heart Fail.* 2014;16:796–801.
62. Nikolic SD, Khairkahan A, Ryu M, et al. Percutaneous implantation of an intra-ventricular device for the treatment of heart failure: experimental results and proof of concept. *J Card Fail.* 2009;15:790–797.